



Most Affordable General Use Seat Cushion



Premier One General Use Foam Cushion

Item #	Product Dimensions	Product HCPCS	UOM
M8061-PB	14" (W) X 14" (D) X 3" (H)	E2601	1/ea, 4/cs
M8062-PB	14" (W) X 16" (D) X 3" (H)	E2601	1/ea, 4/cs
M8063	16" (W) X 16" (D) X 3" (H)	E2601	1/ea, 4/cs
M8064	18" (W) X 16" (D) X 3" (H)	E2601	1/ea, 4/cs
M8065	20" (W) X 16" (D) X 3" (H)	E2601	1/ea, 4/cs
M8070	20" (W) X 18" (D) X 3" (H)	E2601	1/ea, 4/cs
M8077	22" (W) X 18" (D) X 3" (H)	E2602	1/ea, 4/cs
M8078	22" (W) X 20" (D) X 3" (H)	E2602	1/ea, 4/cs
M8079	24" (W) X 18" (D) X 3" (H)	E2602	1/ea, 4/cs
M8080	24" (W) X 20" (D) X 3" (H)	E2602	1/ea, 4/cs

- The Premier One General Use Foam cushion is designed to provide comfort and stability with a simple, economical design.
- Aerated foam helps keep the user cool and dry.
- The low-shear water-resistant stretch cover helps maintain sensitive skin integrity.
- Convenient 12" tie straps secure the cushion to the wheelchair.
- The waterproof vinyl bottom protects the wheelchair seat.
- Meets flammability standard TB117.

specifications	
Cover Material	Fluid-Resistant Stretch Top Waterproof Vinyl Bottom
Weight Capacity	250 lbs. 400 lbs.
Warranty	1 Year Limited